

Department: Radiology	Section: CT	Effective Date: 10/1/2013
Procedure Number: 13.12.16		Revision: 9/2015

Face in Axial Mode 2.9

Clinical Indication: Facial trauma, face pain, cellulitis, tumor

Patient Preparation: None

Orientation: Head first; supine

Breathing: Normal breathing

Oral Contrast: None

IV Contrast Per Weight: *If clinically indicated* 1ml/lb or 2ml/kg Visipaque 320 not to exceed 100ml injected @ 3ml/sec

Coverage: Infraorbital ridge through frontal sinus

Anatomic Reference: Orbitomeatal line

Scan Delay: If indicated use 35 second fixed delay

Group 1:

Scan Mode	Thickness Speed Pitch	Table Interval	SFOV	kV	Auto mA/ Noise Index	Prep Time (sec)	Recon Type
LS 16 1.0 sec Axial Full	5 2i	10	Head	120	280 2.8	N	Standard
VCT 1.0 sec Axial Full	5 4i	20	Head	120	280 2.8	N	Standard
AS 64 2.0 sec Axial	4.8	4.8	200	120	Ref MAS 350	N	Standard

Algorithm: Standard and bone. Recon 1 Standard 2 Bone 3 Standard 0.625mm

Reformat: Sagittal and Coronal

Network: Recon 1 and 2 to PACS; Recon 3 to AWSERV and AWMR