Department: Radiology	Section: CT	Effective Date: 10/1/2013
Procedure Number: 13.12.19		Revision: 9/2015

## Cervical Spine 3.1 Helical

Clinical Indication: Neck pain, fracture

Patient Preparation: None Orientation: Head first Breathing: Normal breathing

**Oral Contrast: None** 

IV Contrast Per Weight: None Coverage: base of skull to T6 Anatomic Reference: EAM

Scan Delay: NA

**Group 1:** If Cervical spine only – base of skull to T6

Scan Mode	Thickness Speed Pitch	Table Interval	SFOV	kV	Auto mA/ Noise Index	Prep Time (sec)	Recon Type
LS 16 1.0 sec Helical Full	2.5 5.62 0.562:1	2.5	Large 16cm	120	<b>380</b> 7	N	Bone Standard
VCT 1.0 sec Helical Full	2.5 5.62 0.562:1	2.5	Large 16cm	120	<b>380</b> 7	N	Bone Standard
AS 64 1.0 sec Helical	2.4	2.4	250	120	Ref MAS 350	N	Bone Standard

Algorithm: Standard and bone. Show recon 1 Standard 2 Bone 3 0.625 @ 25cm Bone

**Reformatting:** Sagittal, coronal and curved coronal off 25cm bone

Network: Recon 1, 2 and reformatted images to PACS; Recon 3 to AWSERV