CT QUALITY CONTROL VISUAL CHECKLIST

GANTRY	Table height indicator								
	Table position indicator								
	Angulation indicator								
	Laser localization inside light								
	Laser localization outside light								
	Table motion smooth								
	X-ray on Indicator								
	Exposure Switch								
	X-ray on indicator								
CONTROL CONSOLE	Display window width/level								
	Panel switches/lights/meters								
	Warning labels								
	Intercom system functioning								
OTHER	Posting "X-ray on light"								
	Posting "Do not enter x-ray in progress								
	Pregnancy notification								
	Service records present								
SMPTE PATTERN	5% patch in 0% square visible								
	95% patch in 100% square visible								
	All greyscale steps visible								
PASS = P	Month								
FAIL = F	Date								
DOES NOT APPLY = NA	Initials								