
 <b>University Health™</b>	<b>POLICY #: Rad Proc 14. 15. 19</b>
<b>SUBJECT:</b> : BODY NON CONTRAST MRA	<b>Effective: 10/1/2013</b> <b>Revised: 2/2015: 10/2017</b>
<b>APPROVED BY</b> <b>Eduardo Gonzalez-Toledo, MD Phd</b>	<b>Page 1 of 2</b>

PLANE	SPC/ LOC	CALIBRA TION STN 3	AX 2D FIESTA ASSET	SAG 2D FIESTA ASSET	SSFSE MPH	INH 3D DELTA UPPER	INH 3D LOWER
SEQ	GRADIE NECHO	GRA ECHO	FIESTA	FIESTA	SPIN ECHO	Inhance Delta Flow	Inhance Delta Flow
MODE	2D	2D	2D	2D	2D	3D	3D
IMAGING OPTIONS	SEQ/ NECHO	FAST/CAL I	SEQ/FAST/ ASSET	SEQ/FAST/ ASSET	Gat/Fast/ MPh/ZIP5 12/SS	FC/Gat/EDR/F AST/IrP/Z512/ ZIP2/ASSET	FC/Gat/EDR/FA ST/IrP/Z512/ZIP 2/ASSET
TE			Min Full	MIN FULL	80	60	60
TR							
FLIP ANGLE			75	75	70		
ETL							70
BW			100	100	83.33	83.33	83.33
FOV	40	48	40	42	40	48	37
SLICE THICKNESS	5	16	7	7	150	2.4	2.4
SLICE SPACING	15	0	2	2	0		
Frequency	256		224	224	320	320	320
Phase	128		320	320	224	224	224
NEX	1		1	1			
PHASE FOV	1		1	1	1	1	1
FREQ DIR	UNSWA P	A/P	R/L	S/I	S/I	S/I	S/I
FLOW COMP DIR							

 <b>University Health™</b>	<b>POLICY #: Rad Proc 14. 15. 19</b>
<b>SUBJECT:</b> : BODY NON CONTRAST MRA	<b>Effective: 10/1/2013</b> <b>Revised: 2/2015: 10/2017</b>
<b>APPROVED BY</b> <b>Eduardo Gonzalez-Toledo, MD Phd</b>	<b>Page 1 of 2</b>

SHIM	AUTO		AUTO	AUTO	AUTO	AUTO	AUTO
PHASE CORRECT	OFF		OFF	OFF	ON	OFF	OFF

ORIENTATION: FEET FIRST/SUPINE      COIL: HD BODYFULL

NOTES: Be sure to include the entire pancreas.

AX SSFSE Thin: Prescribe the slices true axial through the pancreatic duct.

AX LAVA PRE and POST: The mask will be acquired first, then the scanner will pause and the green light will appear. Inject contrast and wait 30 seconds, then start the scan. At 70 seconds, the green light will not appear again. Just click on scan to start this scan. Send the mask, each phase, and the subtractions for each phase to PACS.

OP AX FRFSE Rtr may be utilized if the pt is unable to breathe hold for the exam.