

Note: Additional sequences may be requested at the discretion of the Radiologist monitoring the exam.

EXAM: SHOULDER

ORIENTATION: HEAD FIRST/SUPINE

COIL: SHLDRPA4

PLANE	3 PLN LOC	COR T2	COR T2 FS	SAG T2	SAG STIR	AX GRE	AX T2
SEQ	GRE	FSE XL	FSE XL	FSE XL	FSE IR	GRE	FSE XL
MODE	2D	2D	2D	2D	2D	2D	2D
IMAGING OPTIONS	SEQ/FAST	FC/NPW/EDR/ FAST	FC/NPW/EDR/ FAST	FC/NPW/EDR/ FAST	FC/NPW/SEQ/ EDR/FAST	FC/NPW/EDR	FC/NPW/EDR/ FAST
TE		90 ms +	90 ms +	90 ms +	102 ms	20 ms	90 ms +
TR		3800 ms +	3800 ms +	3800 ms +	3800 ms	500 ms	3800 ms +
TI					150ms		
FLIP ANGLE						20	
ETL		12-18	12-18	12-18	12		12-18
BW		20.83	20.83	20.83	20.83	10.42	20.83
FOV	20	16	16	16	16	16	16
SLICE THICKNESS	7	4	4	4	4	4	4
SLICE SPACING	5	1	1	1	1	1	1
Frequency	256	320	256	320	256	320	320
Phase	128	224	192	224	192	224	224
NEX	1	3	3	3	3	3	3
PHASE FOV	1						
FREQ DIR	UNSWAP	UNSWAP	UNSWAP	UNSWAP	UNSWAP	UNSWAP	UNSWAP
FLOW COMP DIR		FREQ	FREQ	FREQ	FREQ	FREQ	FREQ
SHIM	AUTO	AUTO	AUTO	AUTO	AUTO	AUTO	AUTO
PHASE CORRECT	OFF	ON	ON	ON	ON	ON	ON

Written: 10/2013 Reviewed: 05/2015; 04/2016; 4/2016