

Note: Additional sequences may be requested at the discretion of the Radiologist monitoring the exam.

EXAM: HIP ARTHROGRAM PRE ORIENTATION: FEET FIRST/SUPINE COIL: HD BODY FULL

PLANE	3 PLN LOC	COR STIR	COR T1
SEQ	SPIN ECHO	FSE IR	FSE XL
MODE	2D	2D	2D
IMAGING OPTIONS	SEQ/FAST/SS	FC/NPW/SEQ/EDR/FAST	FC/NPW/EDR/FAST
TE	80	60 ms +	MIN FULL
TR	MIN	3800 ms+	400-600 ms
TI		100-175 ms	
FLIP ANGLE			
ETL		12	2
BW	83.33	20.83	20.83
FOV	44	30	30
SLICE THICKNESS	8	4	4
SLICE SPACING	10	1.5	1.5
Frequency	320	320	320
Phase	192	224	224
NEX		3	3
PHASE FOV	1		
FREQ DIR	UNSWAP	UNSWAP	UNSWAP
FLOW COMP DIR		FREQ	FREQ
SHIM	AUTO	AUTO	AUTO
PHASE CORRECT	OFF	OFF	ON

University Health Shreveport
Department Of Radiology
Proc 14.16.17

Note: Additional sequences may be requested at the discretion of the Radiologist monitoring the exam.

NOTES: Acquire the above sequences and call Dr. Simoncini. Cover the hip of interest.

Written: 05/2014 Reviewed: 02/2015; 04/2016