

Note: Additional sequences may be requested at the discretion of the Radiologist monitoring the exam.

EXAM: HIP ARTHROGRAM POST    ORIENTATION: FEET FIRST/SUPINE    COIL: HD BODY LOWER

PLANE	3 PLN LOC	COR T1 FS C+	AX T1 FS C+	SAG T1 FS C+	AX T1 C+
SEQ	GRE	FSE XL	FSE XL	FSE XL	FSE XL
MODE	2D	2D	2D	2D	2D
IMAGING OPTIONS	SEQ/FAST	FC/NPW/EDR/ FAST	FC/NPW/EDR/ FAST	FC/NPW/EDR/ FAST	FC/NPW/EDR/ FAST
TE		MIN FULL	MIN FULL	MIN FULL	MIN FULL
TR		400-600 ms	400-600 ms	400-600 ms	400-600 ms
TI					
FLIP ANGLE					
ETL		3	3	3	3
BW		20.83	20.83	20.83	20.83
FOV	38	18	18	18	16
SLICE THICKNESS	4	4	4	4	4
SLICE SPACING	5	0.5	0.5	0.5	0.5
Frequency	256	320	320	320	320
Phase	128	224	224	224	224
NEX	1	3	3	3	3
PHASE FOV	1				
FREQ DIR	UNSWAP	S/I	R/L	S/I	R/L
FLOW COMP DIR		FREQ	FREQ	FREQ	FREQ
SHIM	AUTO	AUTO	AUTO	AUTO	AUTO
PHASE CORRECT	OFF	ON	ON	ON	ON

University Health Shreveport  
Department Of Radiology  
Proc 14.16.18

Note: Additional sequences may be requested at the discretion of the Radiologist monitoring the exam.

NOTES: Dr. Simoncini will assist in planning the slices for the post contrast images. Cover the hip of interest.

Written: 05/2014 Reviewed: 02/2015; 04/2016